

Catheter Bag Covers and Tubing Covers

Directions

Catheter Cover

Choose fabric that is machine washable and non-shrinking.

No piecing or design work is necessary. A single fabric is fine to cover or “hide” the catheter bag. As with the pillows we would request male and female friendly themes/fabrics.

The catheter cover is rectangular shaped. It can have rounded or squared corners. The top needs a casing for a “shoestring” type drawstring to be inserted so it can be drawn closed.

Suggested width for rectangular catheter cover is 14 inches wide by 13 inches tall. Two pieces will be needed for front and back or one piece folded to create a 14 x 13 inch bag. The end shape is “purse like”. Three seams are sewn around the edges, with an open casing top for the drawstring.

Tubing Cover

This piece would be separate from the catheter bag but should coordinate with it. The tubing cover should be open at both ends for the tubing to be inserted through it and to be able to attach to the patient on one end and the catheter bag on the other. Please make the tubing cover so that the long side seam is not sewn together but instead “buttoned” together with a series of sew-on Velcro closures. This makes it easy to remove and wash.

Suggested length is 45 inches.

Width of tube is 3 inches

Have 5 – 6 Velcro enclosures along 45 inch seam.

